

## Anti-CCG8 CACNG8 Antibody

Catalog Number: A10890

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




## Overview

|                      |                                                                                                                                            |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-CCG8 CACNG8 Antibody                                                                                                                  |
| Reactive Species     | Human, Mouse, Rat                                                                                                                          |
| Description          | Boster Bio Anti-CCG8 CACNG8 Antibody catalog # A10890. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.      |
| Application          | ELISA, WB                                                                                                                                  |
| Clonality            | Polyclonal                                                                                                                                 |
| Formulation          | Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                        |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit                                                                                                                                     |
| Uniprot ID           | Q8WXS5                                                                                                                                     |

## Technical Details

|                     |                                                                                                                        |
|---------------------|------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Synthesized peptide derived from human protein at AA range: 140-220                                                    |
| Isotype             | IgG                                                                                                                    |
| Form                | Liquid                                                                                                                 |
| Concentration       | 1 mg/ml                                                                                                                |
| Purification        | CCG8 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:500-2000<br>ELISA 1:5000-20000                                                                                    |

## Anti-CCG8 CACNG8 Antibody (A10890) Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4°over night

### Anti-CCG8 CACNG8 Antibody

For Research Use Only. Not for use in diagnostic procedures.