

Anti-IFI44 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A10907-1-Fluoro594

About IFI44

Acts upstream of or within cellular response to virus and response to bacterium. Predicted to be located in cytoplasm. Is expressed in central nervous system and retina. Orthologous to human IFI44 (interferon induced protein 44).

Overview

Product Name	Anti-IFI44 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8BV66

Technical Details

Immunogen	E.coli-derived mouse IFI44 recombinant protein (Position: M1-H370). Mouse IFI44 shares 54.9% amino acid (aa) sequence identity with humanIFI44.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-IFI44 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.