

Anti-GGN Antibody (Center)

Catalog Number: A10914

About GGN

May be involved in spermatogenesis (By similarity).


Overview

Product Name	Anti-GGN Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-GGN Antibody (Center) (Catalog # A10914). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86UU5

Technical Details

Immunogen	This GGN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 139-163 amino acids from the Central region of human GGN.
Predicted Reactive Species	Bovine, Monkey, Mouse
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-GGN Antibody (Center) (A10914) Images



GGN Antibody (Center) western blot analysis in A549, LNCaP cell line lysates (35ug/lane). This demonstrates the GGN antibody detected the GGN protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GGN Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.