

## Anti-IVNS1ABP Antibody Picoband® Cy3 Conjugated

Catalog Number: A10942-2-Cy3

### About IVNS1ABP

Predicted to enable ubiquitin-like ligase-substrate adaptor activity. Involved in RNA splicing; negative regulation of protein ubiquitination; and response to virus. Located in cytosol. Implicated in immunodeficiency 70.

### Overview

Product Name	Anti-IVNS1ABP Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y6Y0

### Technical Details

Immunogen	E.coli-derived human IVNS1ABP recombinant protein (Position: Q30-D580). Human IVNS1ABP shares 97.1% amino acid (aa) sequence identity with mouse IVNS1ABP.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-IVNS1ABP Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.