

## Anti-IVNS1ABP Antibody Picoband® HRP Conjugated

Catalog Number: A10942-2-HRP

### About IVNS1ABP

Predicted to enable ubiquitin-like ligase-substrate adaptor activity. Involved in RNA splicing; negative regulation of protein ubiquitination; and response to virus. Located in cytosol. Implicated in immunodeficiency 70.

### Overview

Product Name	Anti-IVNS1ABP Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y6Y0

### Technical Details

Immunogen	E.coli-derived human IVNS1ABP recombinant protein (Position: Q30-D580). Human IVNS1ABP shares 97.1% amino acid (aa) sequence identity with mouse IVNS1ABP.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IVNS1ABP Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.