

Anti-LUC7L Antibody Picoband® Biotin Conjugated

Catalog Number: A10943-1-Biotin

About LUC7L

Putative RNA-binding protein Luc7-like 1 is a protein that in humans is encoded by the LUC7L gene. The LUC7L gene may represent a mammalian heterochromatic gene, encoding a putative RNA-binding protein similar to the yeast Luc7p subunit of the U1 snRNP splicing complex that is normally required for 5-prime splice site selection.

Overview

Product Name	Anti-LUC7L Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NQ29

Technical Details

Immunogen	E.coli-derived human LUC7L recombinant protein (Position: R78-R362). Human LUC7L shares 96.1% amino acid (aa) sequence identity with mouse LUC7L.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LUC7L Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.