

## Anti-LUC7L Antibody Picoband® Cy3 Conjugated

Catalog Number: A10943-1-Cy3

### About LUC7L

Putative RNA-binding protein Luc7-like 1 is a protein that in humans is encoded by the LUC7L gene. The LUC7L gene may represent a mammalian heterochromatic gene, encoding a putative RNA-binding protein similar to the yeast Luc7p subunit of the U1 snRNP splicing complex that is normally required for 5-prime splice site selection.

### Overview

Product Name	Anti-LUC7L Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NQ29

### Technical Details

Immunogen	E.coli-derived human LUC7L recombinant protein (Position: R78-R362). Human LUC7L shares 96.1% amino acid(aa) sequence identity with mouse LUC7L.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LUC7L Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.