

Anti-KLHL1 Monoclonal Antibody

Catalog Number: A11010

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-KLHL1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-KLHL1 Monoclonal Antibody catalog # A11010. Tested in WB, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 3A8
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9NR64


Technical Details

Immunogen	Purified recombinant fragment of human KLHL1 expressed in E. Coli.
Isotype	IgG1
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-KLHL1 Monoclonal Antibody (A11010) Images



Western Blot analysis using KLHL1 Monoclonal Antibody antibody against KLHL1 recombinant protein.



Immunofluorescence analysis of NIH/3T3 cells using KLHL1 Monoclonal Antibody (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Anti-KLHL1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.