

## Anti-PSTPIP2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A11056-1-Fluoro647

### About PSTPIP2

Proline-serine-threonine phosphatase-interacting protein 2 is an enzyme that in humans is encoded by the PSTPIP2 gene. Predicted to enable actin filament binding activity. Predicted to be involved in actin filament polymerization. Predicted to be located in cytoskeleton and membrane. Predicted to be active in actin filament; cytoplasm; and plasma membrane.

### Overview

Product Name	Anti-PSTPIP2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H939

### Technical Details

Immunogen	E.coli-derived human PSTPIP2 recombinant protein (Position: M1-Q334).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PSTPIP2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.