

Anti-PARVG Antibody Picoband® Biotin Conjugated

Catalog Number: A11061-1-Biotin

About PARVG

Gamma-parvin is a protein that in humans is encoded by the PARVG gene. Members of the parvin family, including PARVG, are actin-binding proteins associated with focal contacts.

Overview

Product Name	Anti-PARVG Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9HBI0

Technical Details

Immunogen	E.coli-derived human PARVG recombinant protein (Position: M1-E280). Human PARVG shares 81.8% amino acid (aa) sequence identity with mouse PARVG.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PARVG Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.