

Anti-Angptl7/Angiopoietin Like 7 Antibody

Catalog Number: A11067-1


Overview

| | |
|----------------------|---|
| Product Name | Anti-Angptl7/Angiopoietin Like 7 Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-Angptl7/Angiopoietin Like 7 Antibody catalog # A11067-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse. |
| Application | ELISA, IF, IHC, WB |
| Clonality | Polyclonal |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | O43827 |

Technical Details

| | |
|---------------------|---|
| Immunogen | The antiserum was produced against synthesized peptide derived from human ANGPTL7. AA range:297-346 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200 |

Anti-Angptl7/Angiopoietin Like 7 Antibody (A11067-1) Images



Western Blot (WB) analysis of SH-SY5Y cells using Angptl7 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Angptl7/Angiopoietin Like 7 Antibody

For Research Use Only. Not for use in diagnostic procedures.