

Anti-Angptl7/Angiopoietin Like 7 Antibody

Catalog Number: A11067-1

About ANGPTL7

Anti-CAND2 antibody is ideal for western blotting, ELISA and IHC. CAND2 is also known as TIP120B, and TATA-binding protein-interacting protein 120B. While both CAND1 (TIP120A) and CAND2 (TIP120B) are TATA-binding proteins and form complexes with various nuclear proteins involved in the control of eukaryotic gene transcription, CAND2 (TIP120B) is expressed specifically in muscle and heart tissue. This is contrary to the ubiquitous expression of CAND1 (TIP120A). TIP120 homologs exist in various higher eukaryotes including *D. melanogaster*, *C. elegans*, and *A. thaliana*. TIP120B is 60% identical in amino acid sequence to TIP120A.

Overview


Product Name	Anti-Angptl7/Angiopoietin Like 7 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Angptl7/Angiopoietin Like 7 Antibody catalog # A11067-1. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal 3A8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43827

Technical Details

Immunogen	Synthesized peptide derived from the C-terminal region of human Angptl7. at AA range: 270-350
Predicted Reactive Species	Bovine, Canine, Guinea Pig, Rabbit
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the

optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
WB, 1:500-1:2000

Anti-Angptl7/Angiopoietin Like 7 Antibody (A11067-1) Images



Western Blot (WB) analysis of SH-SY5Y cells using Angptl7 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Angptl7/Angiopoietin Like 7 Antibody