

Anti-PPEF1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A11074-1-Fluoro594

About PPEF1

Serine/threonine-protein phosphatase with EF-hands 1 is an enzyme that in humans is encoded by the PPEF1 gene. This gene encodes a member of the serine/threonine protein phosphatase with EF-hand motif family. The protein contains a protein phosphatase catalytic domain, and at least two EF-hand calcium-binding motifs in its C terminus. Although its substrate(s) is unknown, the encoded protein has been suggested to play a role in specific sensory neuron function and/or development. This gene shares high sequence similarity with the Drosophila retinal degeneration C (rdgC) gene. Several alternatively spliced transcript variants, each encoding a distinct isoform, have been described.

Overview

| | |
|----------------------|--|
| Product Name | Anti-PPEF1 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | O14829 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human PPEF1 recombinant protein (Position: H137-L652). Human PPEF1 shares 72.7% and 73.2% amino acid (aa) sequence identity with mouse and rat PPEF1, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PPEF1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.