

Anti-POLR1E Antibody Picoband® Biotin Conjugated

Catalog Number: A11077-2-Biotin

About POLR1E

Predicted to enable DNA binding activity and RNA polymerase I general transcription initiation factor binding activity. Involved in nucleolar large rRNA transcription by RNA polymerase I and transcription initiation at RNA polymerase I promoter. Located in fibrillar center and nucleoplasm. Part of RNA polymerase I complex.

Overview

Product Name	Anti-POLR1E Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9GZS1

Technical Details

Immunogen	E.coli-derived human POLR1E recombinant protein (Position: M1-Q291). Human POLR1E shares 78.9% amino acid (aa) sequence identity with mouse POLR1E.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-POLR1E Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.