

Anti-MAST4 Antibody

Catalog Number: A11085-1


Overview

| | |
|----------------------|---|
| Product Name | Anti-MAST4 Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-MAST4 Antibody (Catalog# A11085-1). Tested in WB, IHC-P, IF, ICC, ELISA applications. This antibody reacts with Human, Mouse. |
| Application | ELISA, IF, IHC-P, ICC, WB |
| Clonality | Polyclonal |
| Formulation | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage Instructions | Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | O15021 |


Technical Details

| | |
|---------------------|---|
| Immunogen | A peptide derived from human MAST4. Immunogen sequence location: 251-300 |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Purification | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen. |
| Suggested Dilutions | Western Blot, 1:500 - 1:2000 Immunohistochemistry, 1:100 - 1:300 Immunofluorescence, 1:200 - 1:1000 ELISA, 1:5000 Not yet tested in other applications. |


Anti-MAST4 Antibody (A11085-1) Images



Western blot analysis of lysates from Jurkat, COLO, HUVEC, and MCF-7 cells, using MAST4 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using MAST4 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of COS7 cells, using MAST4 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MAST4 Antibody