

## Anti-GMF Beta/GMFB Antibody Picoband®

Catalog Number: A11087-2

### About GMFB

Predicted to enable Arp2/3 complex binding activity. Predicted to be involved in actin filament debranching and negative regulation of Arp2/3 complex-mediated actin nucleation. Predicted to act upstream of or within learning and locomotory behavior. Predicted to be active in cortical actin cytoskeleton.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-GMF Beta/GMFB Antibody Picoband®                                                                                                                                                                                                                                                                                                                                                                                                         |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Description          | Boster Bio Anti-GMF Beta/GMFB Antibody Picoband® catalog # A11087-2. Tested in WB, ICC/IF, ELISA applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance. |
| Application          | ELISA, IF, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Formulation          | Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .                                                                                                                                                                                                                                                                                                                                                     |
| Storage Instructions | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                             |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Uniprot ID           | P60983                                                                                                                                                                                                                                                                                                                                                                                                                                        |

### Technical Details


|               |                                                                             |
|---------------|-----------------------------------------------------------------------------|
| Immunogen     | E.coli-derived human GMF Beta/GMFB recombinant protein (Position: M1-H142). |
| Form          | Lyophilized                                                                 |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.   |

|                     |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Purification        | Immunogen affinity purified.                                                                                                      |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml, Human, Mouse, Rat<br>Immunocytochemistry/Immunofluorescence, 5 ug/ml, Human<br>ELISA, 0.1-0.5 ug/ml |

## Anti-GMF Beta/GMFB Antibody Picoband® (A11087-2) Images



Western blot analysis of GMF Beta/GMFB using anti-GMF Beta/GMFB antibody (A11087-2). Electrophoresis was performed on a 12% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human U-87MG whole cell lysates, Lane 2: rat brain tissue lysates, Lane 3: mouse brain tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-GMF Beta/GMFB antigen affinity purified polyclonal antibody (A11087-2) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for GMF Beta/GMFB at approximately 17 kDa. The expected band size for GMF Beta/GMFB is at 17 kDa.



IF analysis of GMF Beta/GMFB using anti-GMF Beta/GMFB antibody (A11087-2) and anti-Alpha Tubulin antibody (M03989-3). GMF Beta/GMFB was detected in an immunocytochemical section of U2OS cells. Enzyme antigen retrieval was performed using IHC enzyme antigen retrieval reagent (AR0022) for 15 mins. The cells were blocked with 10% goat serum. And then incubated with 5 ug/mL rabbit anti-GMF Beta/GMFB Antibody (A11087-2) and mouse anti-Alpha Tubulin antibody (M03989-3) overnight at 4°C. Cy3 Conjugated Goat Anti-Rabbit IgG (BA1032) and Fluoro488 Conjugated Goat Anti-Mouse IgG (BA1126) were used as secondary antibody at 1:500 dilution and incubated for 30 minutes at 37°C. Visualize using a fluorescence microscope and filter sets appropriate for the label used.

### Anti-GMF Beta/GMFB Antibody

For Research Use Only. Not for use in diagnostic procedures.