

Anti-PRDM11 Antibody Picoband® Biotin Conjugated

Catalog Number: A11092-2-Biotin

About PRDM11

PRDM11 (PR/SET Domain 11) is a Protein Coding gene. Diseases associated with PRDM11 include Hermansky-Pudlak Syndrome 1. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and methyltransferase activity. An important paralog of this gene is PRDM7.

Overview

Product Name	Anti-PRDM11 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NQV5

Technical Details

Immunogen	E.coli-derived human PRDM11 recombinant protein (Position: E37-L483). Human PRDM11 shares 89.7% amino acid (aa) sequence identity with mouse PRDM11.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PRDM11 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.