

## Anti-PRDM11 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A11092-2-Fluoro647

### About PRDM11

PRDM11 (PR/SET Domain 11) is a Protein Coding gene. Diseases associated with PRDM11 include Hermansky-Pudlak Syndrome 1. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and methyltransferase activity. An important paralog of this gene is PRDM7.

### Overview

Product Name	Anti-PRDM11 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NQV5

### Technical Details

Immunogen	E.coli-derived human PRDM11 recombinant protein (Position: E37-L483). Human PRDM11 shares 89.7% amino acid (aa) sequence identity with mouse PRDM11.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-PRDM11 Antibody - Fluoro647**

**For Research Use Only. Not for use in diagnostic procedures.**