

Anti-ENPP6 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A11099-Fluoro594

About ENPP6

Enables glycerophosphocholine cholinephosphodiesterase activity. Involved in choline metabolic process and lipid metabolic process. Located in extracellular region and plasma membrane.

Overview

Product Name	Anti-ENPP6 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6UWR7

Technical Details

Immunogen	E.coli-derived human ENPP6 recombinant protein (Position: E85-S275).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ENPP6 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.