

Anti-NOL9 Antibody Picoband® Biotin Conjugated

Catalog Number: A11105-1-Biotin

About NOL9

Nucleolar protein 9 is a protein that in humans is encoded by the NOL9 gene. Enables ATP-dependent polydeoxyribonucleotide 5'-hydroxyl-kinase activity. Involved in maturation of 5.8S rRNA. Located in intermediate filament cytoskeleton and nucleolus.

Overview

Product Name	Anti-NOL9 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5SY16

Technical Details

Immunogen	E.coli-derived human NOL9 recombinant protein (Position: E131-N636). Human NOL9 shares 64.8% amino acid (aa) sequence identity with mouse NOL9.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NOL9 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.