

Anti-NOL9 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11105-1-Fluoro550

About NOL9

Nucleolar protein 9 is a protein that in humans is encoded by the NOL9 gene. Enables ATP-dependent polydeoxyribonucleotide 5'-hydroxyl-kinase activity. Involved in maturation of 5.8S rRNA. Located in intermediate filament cytoskeleton and nucleolus.

Overview

Product Name	Anti-NOL9 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5SY16

Technical Details

Immunogen	E.coli-derived human NOL9 recombinant protein (Position: E131-N636). Human NOL9 shares 64.8% amino acid (aa) sequence identity with mouse NOL9.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NOL9 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.