

## Anti-NOL9 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11105-1-Fluoro550

### About NOL9

Nucleolar protein 9 is a protein that in humans is encoded by the NOL9 gene. Enables ATP-dependent polydeoxyribonucleotide 5'-hydroxyl-kinase activity. Involved in maturation of 5.8S rRNA. Located in intermediate filament cytoskeleton and nucleolus.

### Overview

Product Name	Anti-NOL9 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5SY16

### Technical Details

Immunogen	E.coli-derived human NOL9 recombinant protein (Position: E131-N636). Human NOL9 shares 64.8% amino acid (aa) sequence identity with mouse NOL9.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-NOL9 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.