

## Anti-LSM11 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A11132-1-Fluoro594

### About LSM11

LSM11 is a component of the U7 snRNP complex that is involved in the histone 3'-end pre-mRNA processing. Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases.

### Overview

Product Name	Anti-LSM11 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P83369

### Technical Details

Immunogen	E.coli-derived human LSM11 recombinant protein (Position: A17-Q360). Human LSM11 shares 88.4% amino acid (aa) sequence identity with mouse LSM11.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-LSM11 Antibody - Fluoro594**

**For Research Use Only. Not for use in diagnostic procedures.**