

Anti-LSM11 Antibody Picoband® HRP Conjugated

Catalog Number: A11132-1-HRP

About LSM11

LSM11 is a component of the U7 snRNP complex that is involved in the histone 3'-end pre-mRNA processing. Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases.

Overview

Product Name	Anti-LSM11 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P83369

Technical Details

Immunogen	E.coli-derived human LSM11 recombinant protein (Position: A17-Q360). Human LSM11 shares 88.4% amino acid (aa) sequence identity with mouse LSM11.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LSM11 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.