

## Anti-EXOSC7 Antibody Picoband® APC Conjugated

Catalog Number: A11141-1-APC

### About EXOSC7

Predicted to enable mRNA 3'-UTR AU-rich region binding activity. Involved in RNA catabolic process and RNA processing. Located in cytosol and nucleolus. Part of exosome (RNase complex).

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-EXOSC7 Antibody Picoband® APC Conjugated  |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Recommended applications are based on the parent unconjugated antibody (ELISA, IP, IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .   |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.   |
| Host                 | Rabbit   |
| Uniprot ID           | Q15024   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human EXOSC7 recombinant protein (Position: M1-G291).    |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | APC<br>Excitation Wavelength: 633-647 nm<br>Emission Wavelength: 660 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users.                    |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-EXOSC7 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.