

Anti-DCHS2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A11158-Fluoro594

About DCHS2

This gene encodes a large protein that contains many cadherin domains and likely functions in cell adhesion. Genome-wide association studies suggest that this gene may be important in Alzheimer's disease, compressive strength index, and appendicular lean mass.

Overview

Product Name	Anti-DCHS2 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6V1P9

Technical Details

Immunogen	E.coli-derived human DCHS2 recombinant protein (Position: I858-E2983).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-DCHS2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.