

## Anti-SLC7A14 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A11162-2-Fluoro488

### About SLC7A14

This gene is predicted to encode a glycosylated, cationic amino acid transporter protein with 14 transmembrane domains. This gene is primarily expressed in skin fibroblasts, neural tissue, and primary endothelial cells and its protein is predicted to mediate lysosomal uptake of cationic amino acids. Mutations in this gene are associated with autosomal recessive retinitis pigmentosa. In mice, this gene is expressed in the photoreceptor layer of the retina where its expression increases over the course of retinal development and persists in the mature retina.

### Overview

Product Name	Anti-SLC7A14 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8TBB6

### Technical Details

Immunogen	E.coli-derived human SLC7A14 recombinant protein (Position: D9-E771). Human SLC7A14 shares 97% amino acid (aa) sequence identity with mouse SLC7A14.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC7A14 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.