

Anti-KIAA0182/GSE1 Antibody Picoband® HRP Conjugated

Catalog Number: A11186-1-HRP

About GSE1

This gene encodes a proline-rich protein with coiled coil domains that may be a subunit of a BRAF35-HDAC (BHC) histone deacetylase complex. This gene may function as an oncogene in breast cancer and enhanced expression of the encoded protein has been observed in breast cancer patients.

Overview

Product Name	Anti-KIAA0182/GSE1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q14687

Technical Details

Immunogen	E.coli-derived human KIAA0182/GSE1 recombinant protein (Position: L445-A1191).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KIAA0182/GSE1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.