

## Anti-DUSP13/DUSP13B Antibody

Catalog Number: A11187-2

### About DUSP13B

Members of the protein-tyrosine phosphatase superfamily cooperate with protein kinases to regulate cell proliferation and differentiation. This gene encodes a dual specificity phosphatase that acts on both phosphotyrosine and phosphoserine/threonine residues. The encoded protein is expressed in testis.

### Overview

Product Name	Anti-DUSP13/DUSP13B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DUSP13/DUSP13B Antibody catalog # A11187-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UII6

### Technical Details

Immunogen	E.coli-derived human DUSP13/DUSP13B recombinant protein (Position: S3-T195).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DUSP13/DUSP13B Antibody

For Research Use Only. Not for use in diagnostic procedures.