

## Anti-GABRQ Antibody (C-term)

Catalog Number: A11195

### About GABRQ

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.


### Overview

Product Name	Anti-GABRQ Antibody (C-term)
Reactive Species	Human
Description	Boster Bio Anti-GABRQ Antibody (C-term) (Catalog # A11195). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UN88

### Technical Details

Immunogen	This GABRQ antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 585-613 amino acids of human GABRQ.
Predicted Reactive Species	Monkey
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

## Anti-GABRQ Antibody (C-term) (A11195) Images



GABRQ Antibody (C-term) western blot analysis in K562,HL-60,Jurkat,A549 cell line lysates (35ug/lane).This demonstrates the GABRQ antibody detected the GABRQ protein (arrow).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GABRQ Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.