

Anti-Synaptoporin/SYNPR Antibody Picoband® Biotin Conjugated

Catalog Number: A11201-2-Biotin

About SYNPR

Synaptoporin is a protein that in humans is encoded by the SYNPR gene. Predicted to be located in neuron projection and synaptic vesicle. Predicted to be integral component of membrane. Predicted to be active in synaptic vesicle membrane. Predicted to be integral component of synaptic vesicle membrane.

Overview

Product Name	Anti-Synaptoporin/SYNPR Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8TBG9

Technical Details

Immunogen	E.coli-derived human Synaptoporin/SYNPR recombinant protein (Position: M1-I265). Human Synaptoporin/SYNPR shares 95.8% and 98.1% amino acid (aa) sequence identity with mouse and rat Synaptoporin/SYNPR, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Synaptoporin/SYNPR Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.