

## Anti-Synaptoporin/SYNPR Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A11201-2-Fluoro594

### About SYNPR

Synaptoporin is a protein that in humans is encoded by the SYNPR gene. Predicted to be located in neuron projection and synaptic vesicle. Predicted to be integral component of membrane. Predicted to be active in synaptic vesicle membrane. Predicted to be integral component of synaptic vesicle membrane.

### Overview

Product Name	Anti-Synaptoporin/SYNPR Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8TBG9

### Technical Details

Immunogen	E.coli-derived human Synaptoporin/SYNPR recombinant protein (Position: M1-I265). Human Synaptoporin/SYNPR shares 95.8% and 98.1% amino acid (aa) sequence identity with mouse and rat Synaptoporin/SYNPR, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Synaptoporin/SYNPR Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.