

Anti-REM2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A11219-2-Fluoro594

About REM2

Predicted to enable GTP binding activity and calcium channel regulator activity. Predicted to be involved in negative regulation of high voltage-gated calcium channel activity. Predicted to be active in plasma membrane.

Overview

Product Name	Anti-REM2 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8IYK8

Technical Details

Immunogen	E.coli-derived human REM2 recombinant protein (Position: R65-E257). Human REM2 shares 94.8% and 95.9% amino acid (aa) sequence identity with mouse and rat REM2, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-REM2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.