

Anti-STXBP6 Antibody Picoband® HRP Conjugated

Catalog Number: A11292-2-HRP

About STXBP6

Forms non-fusogenic complexes with SNAP25 and STX1A and may thereby modulate the formation of functional SNARE complexes and exocytosis.

Overview

Product Name	Anti-STXBP6 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8NFX7

Technical Details

Immunogen	E.coli-derived human STXBP6 recombinant protein (Position: K9-R176). Human STXBP6 shares 99.4% amino acid (aa) sequence identity with mouse STXBP6.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STXBP6 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.