

## Anti-DNAJC11 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11307-1-Fluoro550

### About DNAJC11

DnaJ homolog subfamily C member 11 is a protein that in humans is encoded by the DNAJC11 gene. DNAJC11 is a member of the J proteins also known as HSP40 family of co-chaperones. Recently DNAJC11 has been found to interact with mitofilin, a mitochondrial inner membrane protein, indicating that it may be a molecular chaperone during mitochondrial protein import and folding.

### Overview

Product Name	Anti-DNAJC11 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NVH1

### Technical Details

Immunogen	E.coli-derived human DNAJC11 recombinant protein (Position: N12-Q551).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNAJC11 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.