

Anti-XTP3TPA/DCTPP1 Antibody Picoband® Biotin Conjugated

Catalog Number: A11334-1-Biotin

About DCTPP1

dCTP pyrophosphatase 1 is an enzyme that in humans is encoded by the DCTPP1 gene. The protein encoded by this gene is dCTP pyrophosphatase, which converts dCTP to dCMP and inorganic pyrophosphate. The encoded protein also displays weak activity against dTTP and dATP, but none against dGTP. This protein may be responsible for eliminating excess dCTP after DNA synthesis and may prevent overmethylation of CpG islands. Three transcript variants, one protein-coding and the other two non-protein coding, have been found for this gene.

Overview

Product Name	Anti-XTP3TPA/DCTPP1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H773

Technical Details

Immunogen	E.coli-derived human XTP3TPA/DCTPP1 recombinant protein (Position: R22-T170).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-XTP3TPA/DCTPP1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.