

Anti-MPP9 MPHOSPH9 Antibody

Catalog Number: A11350-1

About MPHOSPH9

May serve as DNA-binding protein and may be involved in the control of cell-type determination, possibly within the developing nervous system.

Brown L., Proc. Natl. Acad. Sci. U.S.A. 89:8492-8496(1992).

Lipkowitz S., J. Biol. Chem. 267:21065-21071(1992).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Overview

Product Name	Anti-MPP9 MPHOSPH9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MPP9 MPHOSPH9 Antibody catalog # A11350-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99550

Technical Details

Immunogen	Synthesized peptide derived from human TAF5L
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
WB: 1:500-1:1000

Anti-MPP9 MPHOSPH9 Antibody (A11350-1) Images




Figure 1. Western blotting validation for Anti-MPP9 MPHOSPH9 Antibody A11350-1

Western blot (WB) analysis of MPP9 polyclonal antibody at 1:500 dilution

Lane1:Hela cell lysate

Lane2:NIH-3T3 cell lysate

Lane3:PC12 cell lysate

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MPP9 MPHOSPH9 Antibody