

Anti-TFDP3 Antibody (N-term)

Catalog Number: A11396

About TFDP3

This gene encodes a member of the DP family of transcription factors. These factors heterodimerize with E2F proteins to enhance their DNA-binding activity and promote transcription from E2F target genes. This protein functions as a negative regulator and inhibits the DNA binding and transcriptional activities of E2F factors.


Overview

| | |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-TFDP3 Antibody (N-term) |
| Reactive Species | Human |
| Description | Boster Bio Anti-TFDP3 Antibody (N-term) (Catalog # A11396). Tested in WB application(s). This antibody reacts with Human. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q5H9I0 |

Technical Details

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen | This TFDP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 82-108 amino acids from the N-terminal region of human TFDP3. |
| Predicted Reactive Species | Monkey |
| Isotype | Rabbit IgG |
| Purification | This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Suggested Dilutions | WB: 1:1000 |

Anti-TFDP3 Antibody (N-term) (A11396) Images



TFDP3 Antibody (N-term) western blot analysis in WiDr cell line lysates (35ug/lane). This demonstrates the TFDP3 antibody detected the TFDP3 protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TFDP3 Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.