

Anti-IRF2BP1 Antibody Picoband® HRP Conjugated

Catalog Number: A11406-2-HRP

About IRF2BP1

IRF2BP1 (Interferon regulatory factor 2 binding protein 1) is a 64-66 kDa member of the IRF2BP family of proteins. It is widely expressed, and known to bind to the C-terminus of IRF-2, creating a transcriptional repressor complex. IRF2BP1 also bind to JDP2 (Jun dimerization protein 2), promoting its ubiquitination. Human IRF2BP1 is 584 amino acids (aa) in length. It contains two Zn-finger domains (aa 6-67 and 498-584), one coiled-coil region (aa 197-217), a Cys-rich region (aa 503-550), an N-terminal acetylation site at Ala2, and multiple utilized phosphorylation sites at Ser 66/125/186/384/421/436/453/457. Over aa 188-356, human and mouse IRF2BP1 are identical in amino acid sequence.

Overview

Product Name	Anti-IRF2BP1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8IU81

Technical Details

Immunogen	E.coli-derived human IRF2BP1 recombinant protein (Position: A71-D439). Human IRF2BP1 shares 95.1% amino acid (aa) sequence identity with mouse IRF2BP1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-IRF2BP1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.