

## Anti-TRIM45 Antibody

Catalog Number: A11409

### Overview


|                      |  |
|----------------------|--|
| Product Name         | Anti-TRIM45 Antibody   |
| Reactive Species     | Mouse  |
| Description          | Boster Bio Anti-TRIM45 Antibody catalog # A11409. Tested in ICC/IF applications. This antibody reacts with Mouse.                          |
| Application          | IF, ICC  |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q9H8W5   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Recombinant fusion protein of human TRIM45(NP_001139107.1).   |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | ICC/IF: 1:50-1:200  |

## Anti-TRIM45 Antibody (A11409) Images

---



Immunofluorescence analysis of L929 using TRIM45 antibody at dilution of 1:100 . Blue: DAPI for nuclear staining.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-TRIM45 Antibody

For Research Use Only. Not for use in diagnostic procedures.