

Anti-TDRD9 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A11411-Fluoro488

About TDRD9

Predicted to enable ATP hydrolysis activity; RNA binding activity; and helicase activity. Involved in spermatogenesis. Located in cytoplasm and nucleus. Implicated in spermatogenic failure 30.

Overview

Product Name	Anti-TDRD9 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8NDG6

Technical Details

Immunogen	E.coli-derived human TDRD9 recombinant protein (Position: E71-K1339). Human TDRD9 shares 87.6% and 87.3% amino acid (aa) sequence identity with mouse and rat TDRD9, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-TDRD9 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.