

Anti-TDRD9 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11411-Fluoro550

About TDRD9

Predicted to enable ATP hydrolysis activity; RNA binding activity; and helicase activity. Involved in spermatogenesis. Located in cytoplasm and nucleus. Implicated in spermatogenic failure 30.

Overview

| | |
|----------------------|---|
| Product Name | Anti-TDRD9 Antibody Picoband® Fluoro550 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q8NDG6 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human TDRD9 recombinant protein (Position: E71-K1339). Human TDRD9 shares 87.6% and 87.3% amino acid (aa) sequence identity with mouse and rat TDRD9, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-TDRD9 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.