

## Anti-PARP11 Antibody Picoband® APC Conjugated

Catalog Number: A11416-1-APC

### About PARP11

Enables NAD<sup>+</sup> poly-ADP-ribosyltransferase activity and NAD<sup>+</sup>-protein mono-ADP-ribosyltransferase activity. Involved in protein auto-ADP-ribosylation. Located in cytosol; nuclear body; and nuclear envelope.

### Overview

Product Name	Anti-PARP11 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NR21

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human PARP11. Human PARP11 shares 100% amino acid (aa) sequence identity with mouse PARP11.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-PARP11 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.