

Anti-SUPT20H Antibody Picoband® HRP Conjugated

Catalog Number: A11433-2-HRP

About SUPT20H

Predicted to enable transcription coregulator activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within positive regulation of gluconeogenesis and positive regulation of transcription by RNA polymerase II. Part of SAGA-type complex.

Overview

Product Name	Anti-SUPT20H Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8NEM7

Technical Details

Immunogen	E.coli-derived human SUPT20H recombinant protein (Position: M1-Q713).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SUPT20H Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.