

## Anti-DDX19A/B Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A11446-1-Fluoro488

### About DDX19A/B

Predicted to enable RNA helicase activity and mRNA binding activity. Predicted to be involved in poly(A)+ mRNA export from nucleus. Located in membrane.

### Overview

Product Name	Anti-DDX19A/B Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IP, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NUU7/Q9UMR2

### Technical Details

Immunogen	E.coli-derived human DDX19A/B recombinant protein (Position: A9-N478).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-DDX19A/B Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.