

## Anti-TAF11 Antibody

Catalog Number: A11462-1

### About TAF11

TAF11, also named as Transcription initiation factor TFIID 28 kDa subunit, is a 211 amino acid protein, which belongs to the TAF11 family. TAF11 as a core TAFII present in both of the previously described TFIID species which either lack or contain TAFII30. The calculated molecular weight of TAF11 is 23 kDa, but modified TAF11 is about 28kDa. ( PMID: 10336414)

### Overview

Product Name	Anti-TAF11 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TAF11 Antibody catalog # A11462-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q15544

### Technical Details

Immunogen	E.coli-derived human TAF11 recombinant protein (Position: D3-I209).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

---

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TAF11 Antibody

For Research Use Only. Not for use in diagnostic procedures.