

Anti-SH3RF2 Antibody

Catalog Number: A11489-1

About SH3RF2

Enables protein phosphatase 1 binding activity and ubiquitin protein ligase activity. Involved in several processes, including positive regulation of JNK cascade; protein autoubiquitination; and regulation of protein metabolic process. Located in nucleoplasm.

Overview

Product Name	Anti-SH3RF2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SH3RF2 Antibody catalog # A11489-1. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8TEC5

Technical Details

Immunogen	E.coli-derived human SH3RF2 recombinant protein (Position: K77-F292).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SH3RF2 Antibody

For Research Use Only. Not for use in diagnostic procedures.