

Anti-ZNRF2 Antibody

Catalog Number: A11505-2


Overview

Product Name	Anti-ZNRF2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ZNRF2 Antibody catalog # A11505-2. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NHG8

Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 IHC: 1:50-1:200

Anti-ZNRF2 Antibody (A11505-2) Images



Immunohistochemical analysis of paraffin-embedded mouse testis tissue using anti-ZNRF2 antibody. Counter stained with hematoxylin.

1 Publications Citing This Product

1. PubMed ID: -, Qi Xiong,Xiang Tian,Weiling Li,Lin Chen,Mei Zhou,Congyue Xu,Qin Ru Sulforaphane alleviates methamphetamine-induced oxidative damage and apoptosis via the Nrf2-mediated pathway in vitro and in vivo.Food and Agricultural Immunology31(1)2020.https://doi.org

Visit bosterbio.com/anti-znrf2-antibody-a11505-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZNRF2 Antibody

For Research Use Only. Not for use in diagnostic procedures.