

Anti-N-alpha-acetyltransferase 20 NAA20 Antibody

Catalog Number: A11521


Overview

Product Name	Anti-N-alpha-acetyltransferase 20 NAA20 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-N-alpha-acetyltransferase 20 NAA20 Antibody catalog # A11521. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P61599


Technical Details

Immunogen	Recombinant fusion protein of human NAA20(NP_057184.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IP: 1:50-1:200


Anti-N-alpha-acetyltransferase 20 NAA20 Antibody (A11521) Images




Western blot analysis of extracts of various cell lines, using NAA20 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry of paraffin-embedded rat kidney using NAA20 antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded human oophoroma using NAA20 antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunoprecipitation analysis of 200ug extracts of MCF7 cells using 1ug NAA20 antibody. Western blot was performed from the immunoprecipitate using NAA20 antibody at a dilution of 1:1000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-N-alpha-acetyltransferase 20 NAA20 Antibody

For Research Use Only. Not for use in diagnostic procedures.