

Anti-SSH3 Antibody Picoband® Biotin Conjugated

Catalog Number: A11609-3-Biotin

About SSH3

Protein phosphatase Slingshot homolog 3 is an enzyme that in humans is encoded by the SSH3 gene. The ADF (actin-depolymerizing factor)/cofilin family (see MIM 601442) is composed of stimulus-responsive mediators of actin dynamics. ADF/cofilin proteins are inactivated by kinases such as LIM domain kinase-1. The SSH family appears to play a role in actin dynamics by reactivating ADF/cofilin proteins in vivo.

Overview

Product Name	Anti-SSH3 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8TE77

Technical Details

Immunogen	E.coli-derived human SSH3 recombinant protein (Position: Q108-Q620). Human SSH3 shares 82.8% and 84.2% amino acid (aa) sequence identity with mouse and rat SSH3, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SSH3 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.